

# PASS BOX - PASS TROLLEY

## MONOBLOC PASS TROLLEY

## RANGE JG - DECONTAMINATION AIRLOCK



### Description

#### Disinfection pass trolley - Bio-decontamination chamber

Our concept is using a device with high level air-tight that propels the liquid of decontamination as a fog, thus ensuring uniform and rapid action on all the surfaces of the premises to be treated.

The fumigation is an integrated through the wall transfer device for material transfers between different classifications (bio-confinement, ISO classes).

Clear opening : 1500 x 2200 x 2500 mm

Other dimensions on request.

Automatic decontamination with operating screen.

The decontamination cycle is totally automatic with digital signal exchanges between airlock / disinfection system.

Full monitoring of alarms conditions, status functions.

Disinfection operation selected by way of a password protected setup screen on the user screen.

Cycle control and alarm management of pre-defined cycle selectable from stored cycles.

The user interface for starting cycles will be via soft touch buttons on the HMI.

Management: an automaton with touch-sensitive screen allows a use very friendly and offer of multiple functions intended for the use.

#### The management and the safety of functioning.

- Memorization of the made operations.

- Recording of the data – USB key

- Related documentation

- Functional requirement specifications

- Chamber electrical test sign off document

- Chamber commissioning document

- Transfer user manual